

4D-A191 214


AUTOMALYZER SUPPORT FOR USNS LYNCH CRUISE 702-88(U)
TEXAS A AND MUUMIV COLLEGE STATION DEPT OF OCEANOGRAPHY
M A SPEARS ET AL 16 DEC 87 N00014-88-J-6002


1/1

UNCLASSIFIED

F/G 8/3

NL





MICROCOPY RESOLUTION TEST CHART
NATIONAL BUREAU OF STANDARDS-1963-A

AD-A191 214

REPORT DOCUMENTATION PAGE

FILE COPY

1a. REPORT SECURITY CLASSIFICATION UNCLASSIFIED			1b. RESTRICTIVE MARKINGS		
2a. SECURITY CLASSIFICATION AUTHORITY DTIC SELECTED FEB 18 1989			3. DISTRIBUTION AVAILABILITY OF REPORT UNLIMITED		
2b. DECLASSIFICATION/DOWNGRADING SCHEDULE			5. MONITORING ORGANIZATION REPORT NUMBER(S)		
4. PERFORMING ORGANIZATION REPORT NUMBER(S)			7a. NAME OF MONITORING ORGANIZATION Texas A&M Research Foundation		
6a. NAME OF PERFORMING ORGANIZATION Department of Oceanography Texas A&M University		6b. OFFICE SYMBOL (if applicable)	7b. ADDRESS (City, State, and ZIP Code) Box 3578 College Station, Texas 77843		
6c. ADDRESS (City, State, and ZIP Code) College Station, Texas 77843		9. PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER N00014-88-J-6002			
8a. NAME OF FUNDING/SPONSORING ORGANIZATION Naval Ocean R&D Activity (NORDA)		8b. OFFICE SYMBOL (if applicable)	10. SOURCE OF FUNDING NUMBERS		
8c. ADDRESS (City, State, and ZIP Code) Code 334 NSTL, Mississippi 39529		PROGRAM ELEMENT NO	PROJECT NO.	TASK NO	WORK UNIT ACCESSION NO
11. TITLE (Include Security Classification) Autoanalyzer Support for USNS LYNCH Cruise 702-88					
12. PERSONAL AUTHOR(S) M.A. Spears, R.G. Casey & D.C. Biggs					
13a. TYPE OF REPORT TECHNICAL		13b. TIME COVERED FROM 11/1/87 TO 12/31/87		14. DATE OF REPORT (Year, Month, Day) 87-12-16	
15. PAGE COUNT					
16. SUPPLEMENTARY NOTATION					
17. COSATI CODES			18. SUBJECT TERMS (Continue on reverse if necessary and identify by block number)		
FIELD	GROUP	SUB-GROUP	Chemical Oceanography; Nutrient Analyses; Mediterranean Sea; Oceanic Fronts; Chemical Variability in Upper O		
19. ABSTRACT (Continue on reverse if necessary and identify by block number)					
<p>Two Marine Technicians from the Technical Support Services group of the Department of Oceanography at Texas A&M University provided on board autoanalyzer analyses for nitrate, nitrite, phosphate, and silicate on seawater samples collected during Operation DEJA VU, cruise 702-88 of the USNS LYNCH (14 November - 2 December 1987).</p> <p>838 samples from 33 CTD stations were analyzed on board with a Technicon model AA-II four-channel autoanalyzer provided by Texas A&M University.</p> <p><i>George E. Taylor & TAMU</i></p>					
20. DISTRIBUTION/AVAILABILITY OF ABSTRACT <input type="checkbox"/> UNCLASSIFIED/UNLIMITED <input checked="" type="checkbox"/> SAME AS RPT <input type="checkbox"/> DTIC USERS			21. ABSTRACT SECURITY CLASSIFICATION UNCLASSIFIED		
22a. NAME OF RESPONSIBLE INDIVIDUAL Dr. Denis Wiesenburg			22b. TELEPHONE (Include Area Code) 601-688-4600		22c. OFFICE SYMBOL

DISCLAIMER NOTICE

**THIS DOCUMENT IS BEST QUALITY
PRACTICABLE. THE COPY FURNISHED
TO DTIC CONTAINED A SIGNIFICANT
NUMBER OF PAGES WHICH DO NOT
REPRODUCE LEGIBLY.**

SUMMARY OF TECHNICAL SUPPORT SERVICES PROVIDED

Marine Technicians Mark Spears and Glenn Casey of the Department of Oceanography at Texas A&M University flew from Houston, Texas, to Rota, Spain, on 7 November 1987 to provide autoanalyzer support for USNS LYNCH cruise 702-88. Organized by scientists from the Oceanography division of NORDA at Bay St. Louis, Mississippi, this cruise investigated ocean frontal processes in the western Mediterranean Sea, with emphasis on the Almeria and Oran fronts.

Pre-cruise staging and other preparation began in Rota 8 November, and LYNCH embarked Rota on 14 November. From 15 November - 1 December, 41 CTD stations were occupied and four-channel autoanalyzer nutrient analyses were run at 33 of these. Duplicate water samples for analysis of silicate, phosphate, nitrate, and nitrite were drawn at a series of depths from the sample stream from a submersible pump attached to the CTD. Autoanalyzer sampling generally began around 120m and proceeded at roughly 10m intervals to the surface. In all, 419 pairs of water samples were analyzed for nutrients on board LYNCH during cruise 702-88.

SUMMARY OF ANALYTICAL TECHNIQUES

The autoanalyzer used to support LYNCH cruise 702-88 was a 4-channel Technicon AA-II, which was standardized by running 2-3 working standards of all four nutrients prior to and after each station. All samples were analyzed within two hours of collection; analysis rate was 20 samples per hour. Peak heights picked off the strip chart recorder output of each autoanalyzer channel were converted to nutrient concentration in ug-at/liter by linear interpolation from absorbance relative to the working standards, using a shipboard computer.

Silicate was determined by ammonium molybdate + tartaric acid + stannous chloride method; phosphate by ammonium molybdate + hydrazine method; and nitrate by sulfanilamide + NEDA method, after reduction to nitrite with a cadmium reduction column. All analyses except phosphate, which was heated in a 70 C bath, were carried out at room temperature. Colorimeter interference filters utilized were 660 nm (silicate), 880 nm (phosphate) and 550 nm (nitrate + nitrite).

To evaluate precision and accuracy of the analytical methods, 15 replicate analyses of each of the highest working standards were run after completing the last station of the cruise. The standard deviation in ug-at/liter for a 6.0 ug-at/liter silicate standard was 0.05 (cv=1%); 0.005 for a 0.50 ug-at/liter phosphate standard (cv= 1%); 0.05 for a 6.5 ug-at/liter nitrate standard (cv=1%); and 0.005 for a 0.50 ug-at/liter nitrite standard (cv=1%).

TEXAS A & M DEPARTMENT OF OCEANOGRAPHY

PHYSICAL DATA

CRUISE 87L88 STATION B87L88*01*1

WIRE LENGTH	P04	Si(OH)4	N03	N02
METERS	um/l	um/l	um/l	um/l
110	.21	5.80	8.91	100.
100	.21	5.30	6.76	100.
80	.18	3.60	6.60	.02
70	.16	3.00	6.30	.12
60	.12	2.40	5.70	.12
50	.10	1.60	4.13	.19
45	.10	2.20	3.43	.19
40	.05	.50	.55	.03
30	100.	.30	.26	.02
20	.01	.20	.42	.02
10	.02	.10	.65	.01
3	100.	0.00	.40	100.

CRUISE 87L88 STATION B87L88*02*1

WIRE LENGTH	P04	Si(OH)4	N03	N02
METERS	um/l	um/l	um/l	um/l
100	.22	5.10	8.69	.05
90	.22	4.90	9.20	.04
80	.19	4.10	7.24	.05
70	.18	3.80	7.19	.04
60	.17	2.90	5.64	.14
50	.04	.80	4.01	.07
40	.01	.30	.20	.03
30	.01	.20	.20	.01
20	100.	.10	.10	.02
10	.01	.10	.22	.02
0	.03	.20	.77	.03

Accession For

NTIS CRA&I

DTIC TAB

Unannounced

Justification

by
on

Accession Code

by

Specimen

A-1 23

CRUISE 87L88 STATION B87L88*03*1

WIRE LENGTH	P04	Si(OH)4	N03	N02
METERS	UM/1	UM/1	UM/1	UM/1
90	.24	2.50	5.12	.05
80	.21	1.70	4.49	.05
70	.19	1.60	3.93	.07
60	.12	.90	2.58	.03
55	.10	.90	2.08	.05
50	.07	.70	1.89	.05
45	.06	.80	1.63	.03
40	.06	1.00	1.56	.03
32	.04	.70	1.58	.02
20	.02	.30	1.45	.01
10	.02	.30	1.37	.01
0	.02	.10	1.24	100.

CRUISE 87L88 STATION B87L88*04*1

WIRE LENGTH	P04	Si(OH)4	N03	N02
METERS	UM/1	UM/1	UM/1	UM/1
100	.29	3.80	7.24	100.
90	.37	5.00	8.66	100.
80	.35	4.80	8.42	100.
70	.32	4.10	7.55	100.
60	.25	2.70	5.53	100.
50	.20	2.00	4.47	.02
45	.16	3.30	3.41	.16
40	.13	1.20	2.61	.19
30	.01	.10	.12	.14
20	.01	.10	.64	100.
10	.02	.20	.33	100.
0	.01	.20	.30	100.

CRUISE 87L88 STATION 887L88*05*1

WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	UM/I	UM/I	UM/I	UM/I
100	.30	6.00	6.61	100.
80	.22	4.60	7.01	100.
70	.29	4.40	6.40	.12
60	.24	4.00	5.53	.24
55	.16	2.80	4.33	.44
50	.19	2.60	3.80	.43
45	.05	1.10	.48	.17
40	.01	.60	.09	.04
30	.08	.70	.18	.06
20	.02	.80	.27	.07
10	.01	.70	.09	.08
0	.03	.60	.04	.18

CRUISE 87L88 STATION 887L88*06*1

WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	UM/I	UM/I	UM/I	UM/I
100	.35	3.20	8.70	100.
80	.33	2.50	7.68	.01
70	.46	2.20	6.92	.02
60	.32	2.00	6.44	.02
50	.26	1.50	5.68	.09
40	.20	1.20	4.49	.25
35	.20	1.30	4.64	.09
30	.09	.50	1.96	.30
25	100.	99.99	.05	.02
20	100.	99.99	99.99	100.
10	100.	.10	.24	100.
0	.02	.10	.25	.01

CRUISE 87L88 STATION B87L88*07*1

WIRE LENGTH	PO4	SI(OH)4	NO3	NO2
METERS	UM/I	UM/I	UM/I	UM/I
95	.33	3.90	8.36	.01
80	.51	5.10	7.39	.02
70	.26	2.10	5.98	.11
60	.24	2.30	5.47	.15
50	.10	1.10	1.98	.21
40	.06	.80	1.14	.17
35	.03	.60	.72	.12
30	.04	.50	.69	.11
25	.02	.50	.52	.08
20	100.	.40	.21	.03
10	100.	99.99	99.99	100.
0	100.	99.99	99.99	100.

CRUISE 87L88 STATION B87L88*15*1

WIRE LENGTH	PO4	SI(OH)4	NO3	NO2
METERS	UM/I	UM/I	UM/I	UM/I
100	.23	1.50	4.94	.63
90	.21	1.60	4.32	.05
80	.14	1.00	2.74	.05
70	.11	.80	2.10	.07
60	.07	.60	1.00	.37
50	.01	.30	.05	.20
40	.01	.10	99.99	.01
30	100.	.10	99.99	100.
20	100.	.10	99.99	100.
10	100.	.10	.17	100.
0	100.	99.99	.31	100.

CRUISE 87L88 STATION B87L88*16*1

WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
97	.30	3.40	7.16	.03
90	.24	2.60	4.81	.07
80	.27	2.90	5.54	.10
70	.04	.80	.22	.19
60	.01	.60	.01	.06
50	.01	.40	.06	.05
40	.01	.40	.06	.05
30	.02	.50	.06	.06
20	.01	.30	.01	.05
10	100.	99.99	.09	100.
0	.01	.30	.44	100.

CRUISE 87L88 STATION B87L88*17*1

WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
100	.09	2.50	2.90	.04
90	.13	1.50	3.73	.14
85	.14	1.40	3.55	.20
80	.20	2.00	4.74	.26
70	.07	1.20	1.39	.27
60	.05	.90	.78	.50
50	.02	.50	.24	.52
40	.01	.10	.04	.03
30	.01	.10	.04	.02
20	.02	.10	.12	.02
10	.01	.10	.04	.01
0	.01	.10	.05	.01

CRUISE 87L88 STATION B87L88*18*1

WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
100	.14	1.70	4.34	.05
90	.13	1.60	3.72	.08
80	.16	1.60	3.99	.16
70	.13	1.60	2.73	.12
60	.08	1.40	1.32	.28
50	.08	1.20	1.39	.40
40	.04	.80	.60	.51
30	.02	.40	.14	.09
20	.02	.20	.17	.03
15	.01	.30	.14	.04
10	.04	.30	.46	.11
0	.01	99.99	.09	.03

CRUISE 87L88 STATION B87L88*19*1

WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
100	.27	3.50	7.22	.03
90	.27	3.00	7.02	.04
80	.24	2.80	6.53	.06
75	.25	2.80	6.50	.06
70	.25	2.60	6.32	.09
60	.13	1.80	4.77	.25
50	.03	.90	.83	.25
40	.03	.70	.07	.20
35	.02	1.00	.06	.05
30	.02	.60	.04	.02
20	100.	.20	.11	.01
10	.03	.20	.22	.02
0	.01	.30	.07	.02

CRUISE 87L88 STATION B87L88*21*1

WIRE LENGTH	P04	Si(OH)4	N03	N02
METERS	um/l	um/l	um/l	um/l
108	.27	3.30	7.49	.02
90	.27	2.99	7.08	.03
80	.23	2.30	6.61	.05
75	.15	1.50	5.85	.05
70	.09	1.10	2.38	.06
60	.08	1.60	2.06	.06
50	.08	1.80	1.81	.06
40	.07	1.98	1.51	.07
35	.07	1.20	1.55	.29
30	.04	1.30	.66	.22
20	.01	.99	.10	.01
10	.02	.10	.68	.01
0	.03	.10	.26	.03

CRUISE 87L88 STATION B87L88*22*1

WIRE LENGTH	P04	Si(OH)4	N03	N02
METERS	um/l	um/l	um/l	um/l
107	.15	2.80	4.34	.02
90	.16	2.40	4.63	.03
75	.16	2.40	4.35	.03
60	.12	2.00	2.93	.04
50	.18	2.40	4.04	.06
45	.09	1.80	1.87	.10
40	.09	1.80	1.94	.15
35	.05	1.40	.73	.29
30	.02	.70	.57	.11
25	.02	.50	.54	.03
20	.01	.40	.54	.06
10	.01	.20	.20	.02
0	.01	.10	.19	.02

CRUISE 87L88 STATION B87L88x23*1				
WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
108	.22	2.70	6.17	.03
90	.20	2.20	6.02	.03
70	.20	2.30	5.06	.04
60	.17	2.40	3.88	.05
55	.13	1.90	2.70	.07
50	.21	2.20	4.43	.09
40	.18	2.00	3.55	.12
35	.15	1.80	2.88	.25
30	.05	.80	.65	.16
25	100.	.10	.05	.02
20	100.	99.99	.03	.01
10	100.	99.99	.15	.01
0	100.	99.99	.04	.01

CRUISE 87L88 STATION B87L88x24*1				
WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
105	.20	2.70	6.30	.04
95	.18	2.20	5.26	.13
90	.08	1.20	2.35	.35
85	.13	1.20	3.39	.17
80	.16	1.30	3.22	.20
70	.14	1.40	3.06	.34
60	.15	1.60	2.96	.54
50	.21	2.00	4.35	.13
40	.19	1.70	3.86	.12
30	.11	.90	2.05	.15
20	.01	.10	.05	.05
10	100.	.10	.04	.02
0	.01	99.99	99.99	.05

CRUISE 87L88 STATION B87L88*25*1

WIRE LENGTH	P04	Si(OH)4	N03	N02
METERS	um/l	um/l	um/l	um/l
105	.28	3.20	6.65	100.
100	.26	2.60	5.41	.01
90	.23	2.00	4.99	.01
80	.20	1.60	4.04	.02
70	.17	1.30	3.44	.06
60	.11	.70	1.98	.14
50	.02	99.99	.16	.13
45	.02	99.99	.15	.12
40	.02	.20	.16	.12
30	.02	.10	.15	.10
20	.01	.10	.17	.10
10	.01	.10	.06	.10
0	.01	.10	.14	.11

CRUISE 87L88 STATION B87L88*26*1

WIRE LENGTH	P04	Si(OH)4	N03	N02
METERS	um/l	um/l	um/l	um/l
104	.11	2.20	3.74	100.
100	.13	1.70	3.86	.01
90	.31	3.70	7.44	.01
80	.33	3.90	7.45	100.
75	.27	2.90	5.46	.04
70	.32	3.20	6.65	100.
60	.27	2.20	5.21	.01
50	.25	1.90	4.74	100.
40	.14	1.10	2.35	.14
30	.03	.10	.18	.08
20	.01	99.99	.12	.05
10	.02	.10	99.99	.05
0	.03	.10	.06	.06

CRUISE 87L88 STATION B87L88*27*1				
WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
104	.13	2.20	3.87	100.
90	.11	2.30	2.81	.04
80	.09	2.30	2.54	.05
65	.29	3.70	6.62	.07
60	.32	3.50	6.16	.11
50	.29	2.90	5.39	.09
40	.25	2.60	4.54	.12
35	.20	1.90	3.42	.23
30	.08	.80	.82	.25
20	.04	.40	.13	.17
10	.03	.20	99.99	.13
0	.01	.30	99.99	.12

CRUISE 87L88 STATION B87L88*28*1				
WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
109	.27	3.60	6.79	.01
90	.30	3.20	6.73	.01
80	.28	2.80	6.53	.02
70	.27	2.60	6.24	.02
60	.25	2.40	5.32	.02
50	.11	1.20	2.36	.04
40	.24	2.30	4.94	.12
35	.19	2.00	3.60	.12
30	.16	1.90	1.97	.09
25	.01	.30	99.99	.02
20	100.	.10	99.99	100.
10	100.	.10	99.99	100.
0	100.	.10	99.99	100.

CRUISE 87L88 STATION B87L88*29*1

WIRE LENGTH PO4	SI(OH)4	NO3	NO2
METERS	UM/1	UM/1	UM/1
105	.28	2.80	6.67 .03
95	.28	2.70	6.45 .04
90	.29	2.90	6.59 .03
80	.29	2.80	6.54 .04
70	.28	2.70	6.33 .05
60	.25	2.30	5.63 .05
55	.11	1.80	3.60 .07
50	.04	1.00	.44 .12
45	.03	.30	.02 .01
40	.03	.10	.01 .01
35	.02	.20	.02 100.
30	.02	.20	.04 .01
20	.02	.20	.04 .01
10	.03	.20	.04 .01
0	.01	.20	.04 .01

CRUISE 87L88 STATION B87L88*30*1

WIRE LENGTH PO4	SI(OH)4	NO3	NO2
METERS	UM/1	UM/1	UM/1
109	.18	1.80	3.35 .01
100	.33	2.70	5.13 100.
90	.33	2.80	5.03 .01
80	.31	2.20	4.37 .01
70	.29	1.90	3.95 .01
60	.25	1.50	3.22 .02
55	.22	1.30	3.01 .02
50	.21	1.20	2.73 .03
45	.18	1.10	2.38 .03
40	.12	.70	1.55 .12
35	.03	.30	.31 .09
20	.01	.10	.07 .03
10	.01	.20	.07 .03
0	.01	.20	.06 .03

CRUISE 87L88 STATION B87L88*31*1

WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
107	.24	2.10	4.77	.02
100	.24	2.10	4.57	.02
90	.29	2.40	5.02	.02
80	.31	2.50	5.09	.02
70	.30	2.20	4.59	.02
60	.27	2.30	4.21	.02
50	.18	1.40	2.74	.07
45	.19	1.00	1.61	.08
40	.05	.30	.31	.06
35	.05	.20	.10	.04
30	.02	.10	.01	.02
20	.02	.10	.03	.01
10	.02	.10	.02	.01
0	.01	99.99	.02	.01

CRUISE 87L88 STATION B87L88*32*1

WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
107	.10	2.60	5.31	.03
90	.14	2.50	4.33	.03
80	.13	2.50	3.64	.03
70	.08	2.30	2.00	.04
65	.06	1.80	.95	.05
60	.04	1.30	.32	.08
55	.07	1.40	1.37	.15
54	.09	1.40	1.82	.15
50	.03	.80	.11	.10
45	.08	1.40	1.55	.11
40	.03	.70	.15	.08
30	.02	.20	.02	.01
20	.02	.10	99.99	100.
10	.01	.10	99.99	100.
0	100.	99.99	.02	.01

CRUISE 87L88 STATION B87L88*33*1

WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
108	.26	2.70	6.68	.04
90	.24	2.30	6.36	.05
80	.23	2.00	5.66	.06
70	.21	1.60	5.06	.08
60	.23	2.20	5.25	.07
55	.24	2.40	5.14	.06
50	.15	1.70	2.80	.07
45	.13	1.40	2.11	.08
40	.16	1.70	2.58	.09
30	.01	.10	99.99	.01
20	.02	.10	99.99	.01
10	.02	.10	99.99	.01
0	.01	.10	99.99	.01

CRUISE 87L88 STATION B87L88*34*1

WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
105	.18	2.40	5.35	.03
90	.18	2.10	4.76	.03
80	.11	2.40	2.22	.03
70	.09	2.20	1.32	.04
60	.12	1.70	1.53	.08
55	.11	1.40	1.04	.10
50	.10	1.00	.37	.10
45	.03	.10	.07	.03
35	.03	.10	.06	.03
25	.04	.30	.13	.04
20	.05	.20	.09	.03
10	.04	99.99	99.99	.01
0	.02	99.99	99.99	.01

CRUISE 87L88 STATION B87L88*35*1

WIRE LENGTH P04	Si(OH)4	N03	N02
METERS	um/l	um/l	um/l
109	.24	3.30	7.16
100	.26	3.30	7.12
90	.20	3.10	7.06
80	.28	2.90	6.91
70	.26	2.70	6.43
60	.24	2.40	5.78
50	.20	2.00	4.99
40	.14	1.40	3.06
35	.06	.60	.33
30	.04	.20	99.99
25	.03	.40	.01
20	.01	.10	.06
10	.01	.10	.06
0	.01	.10	.06

CRUISE 87L88 STATION B87L88*36*1

WIRE LENGTH P04	Si(OH)4	N03	N02
METERS	um/l	um/l	um/l
107	.11	1.80	2.65
100	.28	3.50	6.21
90	.44	4.60	7.81
80	.45	4.00	6.87
70	.46	3.60	6.16
60	.33	2.40	3.64
50	.10	1.20	.30
45	.08	.80	.16
40	.07	.40	.04
30	.08	.50	.03
20	.05	.50	.02
10	.06	.50	.05
0	.06	.50	.04

CRUISE 87L88 STATION B87L88*37*1

WIRE LENGTH	P04	Si(OH)4	N03	N02
METERS	um/l	um/l	um/l	um/l
107	.14	1.80	3.33	.05
90	.27	2.70	5.97	.04
80	.29	2.90	6.02	.03
70	.24	2.30	4.78	.03
60	.26	2.60	4.72	.04
50	.07	.70	.19	.06
40	.04	.40	.02	.03
30	.03	.30	99.99	.02
20	.02	.30	99.99	.02
10	.02	.30	99.99	.02
0	.02	.30	99.99	.02

CRUISE 87L88 STATION B87L88*36*1

WIRE LENGTH	P04	Si(OH)4	N03	N02
METERS	um/l	um/l	um/l	um/l
108	.23	2.80	6.87	.07
95	.21	2.90	4.59	.06
85	.19	2.30	2.84	.09
80	.15	1.90	1.58	.11
70	.05	.70	99.99	.20
60	100.	.10	99.99	.04
50	100.	.10	99.99	.04
40	.01	.10	99.99	.04
35	.01	.10	99.99	.04
30	.02	.30	99.99	.04
20	.02	.20	.01	.05
10	.02	.20	.01	.04
0	.01	.20	.01	.03

CRUISE 87L88 STATION B87L89*39*1				
WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
109	.14	2.20	4.26	.06
95	.11	2.40	2.01	.05
85	.09	2.60	.74	.04
80	.19	2.40	3.60	.04
70	.12	2.30	.65	.04
60	.13	2.40	.72	.07
55	.02	1.00	.01	.03
50	.01	.30	99.99	.03
40	100.	.30	99.99	.03
30	100.	.30	99.99	.03
20	100.	.20	99.99	.02
10	100.	.10	99.99	.02
0	100.	.30	99.99	.03

CRUISE 87L88 STATION B87L88*40*1				
WIRE LENGTH	PO4	Si(OH)4	NO3	NO2
METERS	um/l	um/l	um/l	um/l
110	.24	3.00	7.31	.08
100	.25	2.90	7.13	.08
90	.24	2.70	6.50	.08
80	.22	2.10	5.27	.07
70	.23	2.30	4.45	.07
65	.15	1.60	1.95	.05
60	.14	1.70	.90	.08
55	.10	1.20	.50	.10
50	.06	.80	.18	.08
40	.04	.40	.01	.07
30	.06	.20	.02	.03
20	.05	.20	.02	.02
10	.03	99.99	99.99	.01
0	.02	.10	99.99	.01

CRUISE 87L88 STATION B87L88*41*1

WIRE LENGTH P04	Si(OH)4	N03	N02
METERS	um/I	um/I	um/I
109	.06	2.60	1.03
95	.08	2.50	1.37
85	.06	2.00	.47
80	.06	1.90	.50
70	.04	.80	.33
60	100.	.10	99.99
50	.01	.10	99.99
40	.02	.10	99.99
30	.02	.10	99.99
20	.01	.10	99.99
10	.01	.10	99.99
0	100.	.10	99.99

DEPTH
IN
METERS

50.0


100

2.00

SILICA UM/1
4.00

6.00

8.00



DEPTH
IN
METERS

50.0


100

2.00

NITRATE $\mu\text{M}/\text{l}$
4.00

6.00

8.00



DEPTH
IN
METERS

50.0

100


PHOSPHATE $\mu\text{M}/\text{l}$

.200

.400


.600

.800



N
I
T
R
A
T
E

U
M
/
l



END

DATE

FILMED

5-88

DTIC